



**7400 Series
Vibrating Sample Magnetometers**





3-year warranty and technical support

Lake Shore products are supported by a 3-year standard warranty, our confirmation of quality and commitment for the long term. Our scientists understand your applications and measurements and provide support throughout your decision making process and beyond the sale.



Vibrating Sample Magnetometers

Contents

System Overview	4 – 5
Materials and Measurements	6
What Our Customers are Saying	7
System Application Software	8 – 9
Post Processing and Analysis	10 – 11
System Pictures	13
7400 Series Specifications	13
7400 Series Equipment	14
Options and Accessories	15 – 17
Options and Accessories	19
Options Specifications	19
Sample Assemblies	20
Sample Rods and Tails	21
Sample Holders Application Notes	22
Shipping Dimensions/Weight	22
Installation Dimensions/Weight	22
Site Requirements	22
Ordering Information	23



Easy sample exchange — the sliding head mechanism allows easy sample exchange and positioning, ensuring reproducibility of measured results


Multiple magnet configurations — 100 mm, 175 mm, and 250 mm (4 in, 7 in, and 10 in) variable-gap electromagnet-based configurations provide fields to 3.1 T

Variable temperature — measure samples from $-269\text{ }^{\circ}\text{C}$ to $1000\text{ }^{\circ}\text{C}$ (4.2 K to 1273 K) with our variable temperature options — the broadest temperature range of any electromagnet-based VSM

Anisotropy measurements — vector coil and autorotation options enable investigations of magnetically anisotropic materials, including derived torque curves

Magneto-resistance probe — perform fast and accurate magneto-resistance measurements with this option as a function of both magnetic field and temperature

Sample holders — accommodate thin film, bulk, liquid, and powder samples





Integrated software — set up and execute measurement routines and experiments quickly and easily from the Windows® menu-driven interface

Detailed post processing — background corrections, automatic offset removal, derivative curves, parameter extraction, and more

Ergonomic workstation — in addition to housing all of the integrated electronics, the workstation acts as a convenient tabletop and has a drawer to store sample holders and samples

Control electronics — the most sensitive electromagnet-based VSM available, featuring a noise floor as low as 1×10^{-7} emu and moment stability of 0.05% per day

The most sensitive, lowest noise floor electromagnet-based VSM system available



Materials

Magnetic thin films, multilayers, and heterostructures

- Particulate
- Continuous
- Magneto-optical
- Magnetic MEMS
- Magnetoresistors (MR)
 - Tunneling-MR (TMR)
 - Giant-MR (GMR)
 - Colossal-MR (CMR)
- Nanomagnetic materials
- Diluted magnetic semiconductors (DMS)
- Paramagnets
- Diamagnets
- Superconductors
- Spin glasses
- Molecular magnetic materials
- Nanocrystalline magnetic alloys
 - Amorphous magnets
 - Melt spun ribbons
- Rare-earth permanent magnets
- Ferrites
 - Hard
 - Semi-hard
- Ferrofluids
- Biological and biomedical
 - Stents
 - MRI contrast agents
 - Nanoscale and microscale particles
- Magnetic powders and inks

Direct and derived measurements as a function of field, temperature, and time

Field-dependent measurements

- Major and minor hysteresis loops
- Saturation magnetization (M_{SAT})
- Remanent magnetization (M_{REM})
- Remanent induction
- B_R
- Coercivity (H_C)
- Intrinsic coercivity (H_{CI})
- Slope at H_C (S^*)
- dM/dH derivative curves
- Differential susceptibility at H_C
- Switching field distribution (SFD)
- Flatness
- Squareness ratio (SQR)
- Initial magnetization curve
- 2nd quadrant demagnetization curves
- Maximum energy
- Product (BH_{MAX})
- DC demagnetization (DCD) remanence
- Isothermal (IRM) remanence
- Permeability curves
- Pinned and free layer parameters
- Exchange field
- Magnetic anisotropy and rotational hysteresis
- Vector (anisotropy) measurements (m_x and m_y)
- Torque curves: $\tau = \mu_0 \mathbf{M} \times \mathbf{H} = -\mu_0 M_y H_x \hat{k}$
- Temperature dependent measurements $M(T)$
 - Curie point
 - Blocking temperature
 - Superconducting transitions and more
- Time dependent measurements $M(t)$
 - Magnetic relaxation
 - Magnetic viscosity

Ideal for the most demanding
magnet characterization
applications



What our customers are saying...

“ We use a Lake Shore VSM to characterize the temperature dependence of the properties of soft magnetic amorphous and nanocrystalline alloys. Current fields of study range from the interaction of superparamagnetic nanoparticles embedded in a matrix to the magnetocaloric effect in amorphous alloys. The high resolution of the system has eliminated our dependence on external SQUID equipment for measurements up to moderate magnetic fields. ”

Dr. Victorino Franco,
Dpto. Física de la Materia Condensada, Universidad de Sevilla, Spain

The Lake Shore VSM purchasing experience was the best ever compared to the other pieces of metallurgical test equipment I have acquired. Lake Shore excelled in technical and sales assistance. I was invited out for training and to witness the QC testing of the equipment before it was shipped, and this was a great help to me. Once I received the equipment, Lake Shore was quick to send someone out for the installation and training. It has always been easy for me to get in touch with the technical representative whenever I had a question.

I have not had a single problem with the power supply, electromagnet, computer, and software. The software and control program seems very logical and is easy to interact with.

I have no regrets about my purchase and am very pleased with the support provided by Lake Shore during the purchasing period and thereafter. They get my highest recommendation.

Jon Stinson,
Boston Scientific Scimed

Magnequench is a premier magnetic material producer that offers products with very tight magnetic property tolerances (some as tight as $\pm 0.6\%$). In order to meet these high quality standards, we use seven Lake Shore VSMs throughout our production processes as well as at our Technical (R&D) center.

Magnequench production VSMs are run continuously (24 hours a day, 7 days a week), and Lake Shore VSMs have proven themselves to be very reliable based on this very demanding environment. Magnequench will most certainly look to Lake Shore when purchasing future VSM systems based on current Lake Shore VSM performance, as well as the wonderful customer support that Lake Shore provides.

Don Kirk Senior Project Engineer,
Magnequench International Inc.

“ Our research focuses on the development of novel EM materials and spintronic devices, which relies heavily on the characterization of magnetic properties of various types of nanostructured materials. The Lake Shore VSM is a workhorse with high sensitivity and rapid measurements that meet our demanding needs. The Lake Shore staff has been extremely helpful in supporting the instrument. What service — I truly appreciate their support of education and research. ”

Dr. John Q. Xiao Professor in Physics,
University of Delaware, Newark

“ The sensitivity of the Lake Shore VSM is as good as specified. We were able to measure ultrathin Co films of 4 Å with the area of $\sim 10 \text{ mm}^2$. We had two publications of Co/Pt multilayers in Physical Review B in the past two years. This instrument facilitates collaborations with my colleagues. I'm particularly satisfied with the low temperature capability. We can measure magnetization down to 8 K with high sensitivity. It replaces our need for SQUID magnetometer for many of our projects. ”

Dr. Fengyuan Yang,
The Ohio State University

Use us as a resource!


Our experts can advise you on the optimal system for your applications. To demonstrate the performance of our VSM and to ensure the proper configuration is selected, we can measure one of your actual samples at no charge to you. Get us involved early and benefit from our many years of experience.

System Application Software

The fully integrated IDEAS VSM software uses an intuitive Windows® interface for system operation, data acquisition, and analysis. Select a default experiment profile or customize your own profile to run a virtually unlimited number of experiments.


Set up and execute measurement routines and experiments quickly and easily. All system parameters and functions are controlled for unattended operation and any number of parameters can be automatically extracted from hysteresis loop data. Real-time field-dependent response tracks field changes for accurate curve shape definitions and parameter extraction.

Substrate corrections and backgrounds can be easily subtracted from measurement data. Calculate and display derivative curves, automatically remove offsets, and determine measurement results. Display real-time feedback of processed data in both tabular and graphical form in CGS or SI units.




Moment vs. Field Setup

Set up a field experiment in either continuous or point-by-point mode.



Ramp Generator

Generate a ramp profile based on field, temperature, or angle setup that is optimized to your material and application



Define Sample Parameters

Define sample parameters before or after recording data.

Application notes

Download free from
www.lakeshore.com or
 request at 614-891-2244.
 Full listing on page 22.

Temperature vs. Moment

Set up temperature-dependent measurement parameters.

Put our IDEAS™ application software in control

Specify Adjustments and Calculations

Correct for offsets, sample holder or substrate contributions, and linear background slopes. Calculate derivative curves and determine exchange and pinning fields.


Profile Manager

Define, save, and edit individual experiments as well as versatile profiles. Profiles allow you to automate sequences of multiple experiments along with other parameters, such as start date and time, and field and temperature ramps.




Measurements


11 μemu CoPt Hard Disk Film



Magnetic Tape — Minor Loop Results




Vector Results



graph represents multiple data sets overlaid


Take even the most complex measurements with ease

CoPt Thin Film — $M(H)$ for H parallel and perpendicular to film plane




graph represents multiple data sets overlaid

NdFeB — Initial magnetization, minor and major loops






NiFe thin film (3 nm) — 3.5 nAm² (3.5 μ emu) and 4 A/m (0.05 Oe) Steps




CMR Film — Low Temperature Results




graph represents multiple data sets overlaid

m_y vs. H_x vs. θ




graph represents multiple data sets overlaid

M(T) on warming & cooling for a nanocrystalline melt-spun ribbon (NdFeB_x)




graph represents multiple data sets overlaid

Magnetic Tunnel Junction with Al₂O₃ Barrier*



graph represents multiple data sets overlaid

Magnetic Tunnel Junction with MgO Barrier*



*Thanks to Professor Dr. John Q. Xiao's group at the University of Delaware, Newark for providing the samples for these measurements

**Model 7404****Model 7407****Model 7410**



7400 Series VSM Specifications



	7404	7407	7410
Magnet pole cap face diameter	51 mm (2.0 in) ¹	51 mm (2.0 in) ¹	50 mm (2.0 in) ¹
APPLIED FIELD STRENGTH (±1%)			
Room temperature			
16.2 mm (0.64 in) magnet air gap/ 3.6 mm (0.14 in) sample access	21.7 kOe (2.17 T)	26.2 kOe (2.62 T)	31 kOe (3.1 T)
23 mm (0.9 in) magnet air gap/ 10 mm (0.4 in) sample access	18.1 kOe (1.81 T)	23.4 kOe (2.34 T)	28.5 kOe (2.85 T)
29 mm (1.14 in) magnet air gap/ 16 mm (0.64 in) sample access	15.3 kOe (1.53 T)	21.2 kOe (2.12 T)	25 kOe (2.5 T)
With oven or cryostat option			
38.1 mm (1.5 in) magnet air gap/ 6.4 mm (0.25 in) sample access	11.8 kOe (1.18 T)	18.3 kOe (1.83 T)	23 kOe (2.3 T)
With single-stage variable temperature option			
33 mm (1.3 in) magnet air gap/ 6.4 mm (0.25 in) sample access	13.5 kOe (1.35 T)	19.4 kOe (1.94 T)	25.0 kOe (2.5 T)
MOMENT MEASUREMENT			
Noise floor (emu RMS)			
Room temperature			
0.1 TC; 0.1 s/pt (no averaging)		0.75 μemu	
0.1 TC; 1 s/pt		0.4 μemu	
0.1 TC; 10 s/pt		0.1 μemu	
With oven or cryostat option, 0.1 TC, 10 s/pt		2.5 μemu	
With single stage variable temperature option 0.1 TC, 10 s/pt		1.25 μemu	
With vector option 0.1 TC, 10 s/pt	5 μemu		3 μemu
Dynamic range		1 × 10 ⁻⁷ to 10 ³ emu	
Time constants (TC)		0.1, 0.3, 1.0, 3.0, or 10.0 s	
Moment stability ²	Less than 0.05% RMS of full scale/day at constant field, constant temperature, and highest range after a 5 h warm-up period		
Reproducibility	Better than ±0.5%, or ±0.1% of full scale, fixed rotation angle and range, with sample replacement		
Moment accuracy	Better than 1% of reading ±0.2% of full scale with a geometrically identical test sample and calibrant		
Sample mass	0 to 10 g (higher mass can be accommodated with decreased performance)		
FIELD MEASUREMENT			
Field accuracy	1% of reading or ±0.05% of full scale		
Field resolution			
2800 kA/m (35 kOe)		8 A/m (0.1 Oe)	
280 kA/m (3.5 kOe)		0.8 A/m (0.01 Oe)	
28 kA/m (350 Oe)		0.08 A/m (0.001 Oe)	
Closed loop field control stability	<0.05% RMS of full scale/h		
MANUAL ROTATION			
Setting resolution	<1°		
Setting reproducibility	<1°		
Rotation range	0 to 730°		
CERTIFICATIONS			
CE	yes		
Application of Council directives	73/23/EEC; 89/336/EEC		
Standard to which conformity is declared	EN61010-1: Overvoltage II, Pollution Degree II; EN61326: Class A, Annex B		
UTILITIES			
Total system cooling water power dissipation (50 or 60 Hz) — contact Lake Shore for the most current selection of available recirculating chillers	4.25 kW; 7.45 kW with optional 74018	13.4 kW; 16.6 kW with optional 74018	13.4 kW; 16.6 kW with optional 74018

¹ With standard 740EMSC coils

² Tested with an AlNiCo sample in a 1 inch sensing coil gap after system warm-up period with the sample vibrating at field. The AlNiCo samples' moment must be >50% of full scale moment range.



7400 Series VSM Equipment

	7404	7407	7410
VSM head drive		Model 74014	
VSM frame		1	
Control electronics		Model 736	
Linear amplifier		Model 142	
Bipolar magnet power supply	Model 643		Model 648
Mode		DC current source	
Maximum output	±35 V/±70 A (2450 W)		±75 V/±135 A (9.1 kW nominal)
AC line input	204/8 VAC ±10%, 13 A/phase; 220/230 VAC ±10%, 12 A/phase; 380 VAC ±10%, 7 A/phase; 400/415 VAC ±10%, 6.5 A/phase at 50/60 Hz		200 VAC ±10%, 41 A/phase; 208 VAC ±10%, 40 A/phase; 220 VAC ±10%, 38 A/phase; 230 VAC ±10%, 37 A/phase; 380 VAC ±10%, 23 A/phase; 400 VAC ±10%, 21 A/phase; 415 VAC ±10%, 21 A/phase
Cooling water requirements	Tap water or closed cooling system (optional chiller available) +15 °C to +30 °C		Tap water or closed cooling system (optional chiller available) +15 °C to +30 °C
Flow rate	5.7 L (1.5 gal)/min minimum		7.6 L (2.0 gal)/min minimum
Pressure drop	10 kPa (1.5 psi) at 5.7 L (1.5 gal)/min minimum for power supply only		159 kPa (23 psi) at 7.6 L (2.0 gal)/min minimum for power supply and mandatory flow switch
Electromagnet	Model EM4-HVA	Model EM7-HV	Model EM10-HV
Pole diameter	100 mm (4 in)	180 mm (7 in)	250 mm (10 in)
Pole cap face diameter	50 mm (2 in)	50 mm (2 in)	50 mm (2 in)
Field homogeneity	±0.1% over 1 cm ² (0.4 in ²)	±0.1% over centered 5 cm (2 in) diameter circle	±0.1% over centered 5 cm (2 in) diameter circle
Cooling water requirements	Tap water or closed cooling system (optional chiller available)	Tap water or closed cooling system (optional chiller available)	Tap water or closed cooling system (optional chiller available)
Inlet temperature	15 – 25 °C (59 – 77 °F)	15 – 32 °C (59 – 89 °F)	15 – 25 °C (59 – 77 °F)
Flow rate	7.6 L (2 gal)/min	11.4 L (3 gal)/min	15 L (4 gal)/min
Pressure drop	200 kPa (30 psi)	220 kPa (32 psi)	200 kPa (30 psi)
Water chiller capacity	2.5 kW (8,530 BTU)/h	5 kW (17,060 BTU)/h	8.8 kW (30,035 BTU)/h
Hall probe	High stability; 74 mm (2.9 in) aluminum stem		High stability; 203 mm (8 in) aluminum stem
Instrument console		483 mm (19 in) rack	
Computer with IDEAS™ software		Model 740935	

The chillers we offer are rated at 65% duty cycle. This is appropriate for many common magnet testing applications such as hysteresis loops and other measurement applications where the operating cycle is spent at low to medium current, with only limited excursions to high fields. It may be appropriate to choose a larger chiller for higher duty cycle needs, when large magnetic fields must be maintained continuously. Please consult Lake Shore for an appropriate chiller for these applications.



Expand your capabilities with options

Model 74035 single-stage variable temperature option

The single stage variable temperature assembly allows you to take measurements from 100 K to 950 K using LN₂, nitrogen, and argon gas. A single point measurement can be taken at 78 K. Only one hardware device is required to go from high to low temperatures, eliminating the need to remove or resaddle your sample. This ensures accurate measurements throughout the full scale temperature range. Rapid cool down from 950 K to room temperature and from room temperature to 100 K provides efficiency and high throughput. Like our full suite of variable temperature options, the single stage variable temperature assembly is mechanically isolated from the magnetometer head and sample, minimizing noise floor. Designed to deliver superior thermal performance, the unit's vacuum insulation prevents freeze over at low temperatures and can operate safely at high temperatures without the risk of damaging neighboring components.

Included with Model 74035:

1. Single-stage variable temperature insert with mount
2. 25 liter LN₂ Dewar with condenser stand
3. Gas handling box
4. LN₂ transfer line with condenser assembly
5. Instrument cables
6. Sample rods and holders
 - a. With 7404 and 7407: 740928 sample tail and holder kit
 - b. With 7410: 740941 sample tail and holder kit

Supplemental 74035 equipment requirements:

1. Lake Shore Model 741-VTA variable temperature option kit
2. Argon gas cylinder with 344 kPa (50 psi) gas regulator and 3 mm hose barb (can also be a ¼ NPT female fitting)
3. Nitrogen gas cylinder with a 344 kPa (50 psi) gas regulator and 3 mm hose barb (can also be a ¼ NPT female fitting)
4. LN₂ source to fill the provided Dewar
5. Clean compressed air (276 kPa [40 psi])
6. Mechanical vacuum pump (E2M or equivalent) kit providing sample space blank off pressure of 0.67 Pa (5×10^{-3} torr) for routine operation
7. Turbomolecular vacuum pump (Lake Shore TPS-FRG or equivalent) kit for cryogen transfer line maintenance — can also be used in place of the E2M rotary vacuum pump
8. A Pirani or thermocouple vacuum gauge capable of measuring pressures from 0.1 to 100 Pa (10^{-3} to 1 torr)

Model 74018 variable temperature cryostat

The 7400 Series VSM cryostat is designed for rapid sample cooling with either LHe or LN₂ as well as easy sample insertion and interchange. It allows you to take measurements from 5.5 K to 450 K using LHe and from 85 K to 450 K using LN₂. A single-point measurement can be taken at 4.2 K (LHe) and at 77.6 K (LN₂).

The sample is suspended in a proprietary insulated tube constructed of nonmagnetic material.

The cryostat is mechanically isolated from the magnetometer head and sample, greatly reducing the system noise floor.

It is mounted between an electromagnet base plate and a quick release mechanism located on the top of the electromagnet.

The cryostat design provides the user the capability to perform measurements economically over nearly the entire accessible temperature range with a single cryostat. The transfer line is included with the cryostat.

Included with Model 74018:

1. Combination LHe/LN₂ cryostat with mount
2. LHe/LN₂ transfer line
3. Cryogen transfer kit
4. Instrument cables and related accessories
5. Sample rods and holders
 - a. With 7404 and 7407: 740929 sample tail and holder kit
 - b. With 7410: 740943 sample tail and holder kit

Supplemental 74018 equipment requirements:

1. Lake Shore Model 741-VTA temperature option kit
2. A mechanical vacuum pump (Lake Shore E2M or similar) capable of achieving a pressure below 0.67 Pa (5×10^{-3} torr) and a speed of 1 m³/h, along with a KF-16 flange pump inlet
3. Access to turbomolecular vacuum pump (Lake Shore TPS-FRG or similar) capable of doing better than 1.33×10^{-3} Pa (10^{-6} torr) for annual evacuation of transfer line vacuum space
4. LHe or LN₂ storage Dewar (Lake Shore 1220-50 or similar) with top withdraw fitting to accept the 12.7 mm (0.5 in) diameter transfer line — the transfer line furnished with the Model 74018 cryostat is particularly well adapted for use with 25 to 60 L storage vessels, and can be readily adapted to other capacity storage vessels (in most cases, a LHe Dewar will be provided by your local liquid gas distributor when LHe is delivered)
5. Gas cylinder with 1 to 5 psi pressure regulator to deliver clean, dry helium or nitrogen gas (depending on liquid cryogen)



Configure your system —

Model 74034 high temperature oven

The high temperature oven assembly enables the Model 7400 VSM system to be used to investigate the magnetic properties of materials at high temperature. This option consists of an electrically heated outer tube assembly with efficient thermal insulation to permit sample-zone temperature from 100 °C to 1000 °C (373 K to 1273 K). Temperatures from 30 °C to 1000 °C (303 K to 1273 K) are also possible, however, below 100 °C (373 K) measurement time increases.

The inner sample zone chamber is lined with a special heat-resistant and intrinsically non-magnetic material. A sample holder is provided which consists of a quartz tube sample rod attached to a boron-nitride sample cup. The oven secures to a special isolation mount support structure, the only special consideration being that the air gap between the coils must accommodate the 23 mm (0.9 in) outside diameter of the oven. A mechanical vacuum pump capable of maintaining inlet pressures down to 0.67 Pa (5×10^{-3} torr) must be supplied by the user.

This option features efficient thermal insulation, consisting of an evacuation outer chamber with multiple reflective heat shields. Sample zone temperatures as high as 1000 °C are attained with a power consumption of approximately 70 W. Two results of the low power consumption are minimal magnetic interference and increased temperature uniformity in the sample zone. The oven is particularly well suited to measuring Curie temperatures of ferromagnetic or ferrimagnetic materials at temperatures up to 1000 °C. The sensitivity of the Model 7400 VSM permits Curie temperature determinations at relatively low field intensities, allowing more inherently accurate determinations.

At room temperature and above, measurements may be performed on samples contained in an air or argon atmosphere to protect the sample from oxidation.

Included with Model 74034:

1. Oven assembly with mount
2. Gas handling box
3. Nickel Curie sample cylinder
4. Instrument cables and related accessories
5. Sample rods and holders
 - a. With 7404 and 7407: 740928 sample tail and holder kit
 - b. With 7410: 740941 sample tail and holder kit

Supplemental 74034 equipment requirements:

1. Lake Shore Model 741-VTA temperature option kit
2. A mechanical vacuum pump (Lake Shore E2M or similar) capable of achieving a blanked-off pressure below 0.67 Pa (5×10^{-3} torr) and a pumping speed of 1 m³/h, along with a KF-16 flange pump inlet
3. Argon gas cylinder with 5 to 10 psi regulator and 3 mm (1/8 in) hose barb



Model 741-VTA temperature option kit

The autotuning cryogenic temperature controller is used to measure and control our full suite of variable temperature options. The Model 741-VTA includes a Lake Shore temperature controller, thermocouple



input card (when purchased for use with the high temperature oven or single stage variable temperature assembly), vacuum handling kit, mounting hardware (included with Model 74035), flanges, hoses, connectors, and accessories. Note: only one 741-VTA is required for all variable temperature options.



exactly the way you need it

Model 74046 magnetoresistance (MR) probe

The MR probe option performs fast and accurate measurements of MR, GMR spin-valve, CMR and other magnetoresistive materials as a function of both in-plane magnetic field and temperature. This measurement option includes data acquisition, control, and analysis software to automatically extract pertinent parameters for the device under test. These include free and pinned layer parameters of both simple and synthetic spin-valve sensors.

The MR probe features four in-line pins for solder-less connection to the sample and is interchangeable with the VSM sample rod for compatibility with variable temperature and autorotation options. The temperature range of use extends from 20 K to 450 K when used with the 74018 LHe/LN₂ cryostat, 325 K to 673 K with the 74033 oven, and 100 K to 673 K with the 74035 single-stage variable temperature assembly. Angular dependent MR measurements are possible when used in combination with the 74033 autorotation option. Contact pins maintain their position on the sample while the entire assembly rotates with respect to the magnetic field.

The fully automated MR software is an intuitive yet powerful user interface providing automatic control of all experimental parameters for unattended operation. Experiment recipes can be saved, retrieved, and edited, and measurement data can be displayed and exported in graphical or tabular format. Multiple step profiles can also be defined to allow for flexibility in the preparation steps and for developing annealing step process methodologies.

The MR probe option is composed of several user-replaceable parts, including contact pins, ceramic pin guides, a ceramic sample holder, a contact pressure spring, and a printed circuit board, providing the convenience of in-field maintenance.

Included with Model 74046:

1. MR hardware insert
2. Lake Shore Model 776 matrix switch
3. Model 2400 Keithley sourcemeter
4. 4-wire I-V cable and MR adapter box
5. MR spare kit (2 PCBs, 1 contact pressure spring, 8 points, 1 pin guide, 1 lower holder — user-replaceable)



Model 74033 autorotation

The rotation option allows you to automatically vary the sample orientation relative to the direction of the applied magnetic field. The angle of rotation is within a single plane defined by the direction of applied magnetic field, referred to as the x-axis. Angular variation is about the z-axis. Rotation is programmable to a resolution of <math><1^\circ</math> for rotating the sample from -10 to 730° and all parameters are measured as a function of rotation angle.

Model 74032 vector coils

The vector option extends the VSM measurement capabilities to facilitate investigations of anisotropic magnetic materials, allowing you to determine their vector magnetization components and susceptibility tensor. When used in combination with the Model 74033 autorotation option, the vector coils provide information that is essentially identical to that provided by a dedicated torque magnetometer. The Model 74032 2-inch vector coils are compatible with all variable temperature options.

E2M 2-stage rotary vacuum pump

Capable of achieving a pressure below 0.67 Pa (5×10^{-3} torr) at 1 m³ per hour, the 2-stage rotary vacuum pump is used for evacuating both the oven and cryostat vacuum spaces of our variable temperature options. This, or a similar vacuum pump, is required for daily operation of variable temperature options.

TPS-FRG turbomolecular vacuum pump station

Used to annually evacuate the cryogen transfer line of the optional cryostat and single stage variable temperature assembly (transfer line and kit are included with these options), the Model TPS-FRG provides vacuum to 1.33×10^{-3} Pa (10^{-6} torr). In addition to annual cryogen transfer line maintenance, the turbomolecular vacuum pump can also be used in place of the E2M rotary vacuum pump for evacuating the cryostat vacuum space.

Recirculating chillers

Lake Shore offers NesLab[®] recirculating chillers in order to provide a complete laboratory solution. The NesLab chillers feature a CFC-free refrigeration system.

The refrigeration system utilizes a hermetically sealed compressor and hot gas bypass system of temperature control. This system eliminates on/off cycling and premature wear of the compressor. Strong pumps provide continuous flow even through cooling lines with small IDs.



7400 Series VSM Option Specifications

Model 74018 variable temperature cryostat

With LHe	Temp range	4.2 K base, 5.5 K to 450 K control
	Temp stability	±0.1 K
With LN ₂	Temp range	77.6 K, 85 K to 450 K
	Temp stability	±0.2 K
Temperature resolution		0.001 K
Cool-down time		5 min (15 min initial cool-down)
Nominal ramp rate		Continuous flow
LHe liquid usage		<1 L/h when operating >7 K
LN ₂ liquid usage		<1 L/h when operating >7 K
Insulation		Vacuum
Sample zone dimensions	Bore size	7.1 mm (0.28 in)
	Outside diameter	22.4 mm (0.88 in)

Model 74033 autorotation

Full range of rotation	-10 to 730°
Setting resolution	<1°
Setting repeatability	<1°

Model 74034 high temperature oven

Temp range	303 K to 1273 K	
Temp stability	±0.1 K	
Temp resolution	0.001 K	
Nominal ramp rate	353 K/min at maximum heating rate of 80 W	
Insulation	Vacuum plus multiple reflective shields	
Sample zone dimensions	Bore size	7.1 mm (0.28 in)
	Outside diameter	25 mm (1.0 in)

Model 74035 single stage variable temperature option

Temperature range	78 K (base), 100 K to 950 K (control)	
Temperature stability	±0.1 K	
Temperature resolution	0.001 K	
Gasses	LN ₂ and nitrogen gas for T<350 K; argon for T>350 K	
Cool-down time	15 min from room temp to 100 K, 40 min from 1000 K to room temp	
Nominal ramp rate (in the domain)	5 K/min	
Hold time	Continuous flow	
LN ₂ usage	<0.75 L/h >100 K – 350 K	
Nitrogen gas usage	3.2 L/min 100 K – 350 K	
Argon gas usage	3.6 L/min	
Insulation	Vacuum	
Sample zone dimensions	Bore size	7.1 mm (0.28 in)
	Outside diameter	17.8 mm (0.7 in)

Model 74033 autorotation

Full range of rotation	-10 to 730°
Setting resolution	<1°
Setting repeatability	<1°

Model 74046 magnetoresistance (MR) probe

Number of probes	4
Probe pin to pin spacing	0.9 mm
Total 4-pin spacing	2.7 mm
Nominal sample size	4.5 mm × 4.5 mm cross-section, maximum height 3 mm
Temperature range	20 K to 673 K
Current ranges	6 ranges; 1 μA to 100 mA
Resistance ranges	9 ranges: 0.2 Ω, 2 Ω, 20 Ω, 200 Ω, 2 kΩ, 20 kΩ, 200 kΩ, 2 MΩ, up to 10 MΩ
Probe tip compliance voltage	0 V to 5 V, measurement 0 V to 100 V, contact formation

	7404	7407	7410
Model 74032 standard vector coils⁴			
Pole caps	102 mm (4 in) ⁵	102 mm (4 in) ⁵	100 mm (4 in) ⁵
Air gap	50.8 mm (2 in)	50.8 mm (2 in)	50.8 mm (2 in)
Maximum applied field	7.7 kOe (0.77 T)	12.5 kOe (1.25 T)	20 kOe (2.0 T)
RMS noise (Y coils)	5 μemu	5 μemu	3 μemu
Sample access		25 mm (1 in)	
Torque density minimum (at maximum applied field)	38 × 10 ⁻³ dyn-cm	62 × 10 ⁻³ dyn-cm	60 × 10 ⁻³ dyn-cm

⁴ With 730ESC coils; purchased separately

⁵ Pole caps sold separately



Sample Assemblies

Lake Shore provides sample rod assemblies for room, cryogenic, and high temperature applications. The assemblies are available as permanently mounted integrated sample tail/holders or with threaded ends that allow holders to be interchanged. Sample holders for thin films, bulk materials, powders, and liquids are available.

The following charts show which sample kits are included standard with different VSM system configurations.

7404 and 7407 System Sample Kits

SYSTEM	KIT	Integrated Sample Tail/Holders				Sample Tails		Sample Holders						
		740930	740931	740933	740934	740932	740935	730931	730933	730934	730935	730937	730938	730939
		Quartz, thin film bottom	Quartz, thin film side	Fiberglass, 3.5 mm air gap, thin film side	Fiberglass, 3.5 mm air gap, thin film bottom	Quartz to BN	Fiberglass	Kel-F bulk/powder upper and bottom cup	Kel-F thin film side	Kel-F thin film bottom	Kel-F liquid upper and bottom cup	Disposable BN cup	BN thin film side	BN thin film bottom
Room temperature	740927			1	1			1	3	3	3	1		
With 74034 variable temperature oven	740928	1	1			1						1	1	1
With 74018 variable temperature cryostat	740929	1	1											
With 74035 single stage variable temperature	740928	1	1			1						1	1	1


7410 System Sample Kits

SYSTEM	KIT	Integrated Sample Tail/Holders				Sample Tails		Sample Holders						
		740937	740938	740944	740945	740939	740942	730931	730933	730934	730935	730937	730938	730939
		Quartz, 3.5 mm air gap, thin film side	Quartz, 3.5 mm air gap, thin film bottom	Quartz, thin film bottom	Quartz, thin film side	Quartz to Kel-F	Quartz to BN	Kel-F bulk/powder upper and bottom cup	Kel-F thin film side	Kel-F thin film bottom	Kel-F liquid upper and bottom cup	Disposable BN cup	BN thin film side	BN thin film bottom
Room temperature	740940	1	1			1		3	3	3	1			
With 74034 variable temperature oven	740941			1	1		1					1	1	1
With 74018 variable temperature cryostat	740943			1	1									
With 74035 single stage variable temperature	740941			1	1		1					1	1	1


Integrated Sample Tail/holders

For Model 7404/7407

740930 Quartz, thin film bottom



740931 Quartz, thin film side



740933 Fiberglass, 3.5 mm air gap, thin film side




740934 Fiberglass, 3.5 mm air gap, thin film bottom




For Model 7410


740937 Quartz, 3.5 mm air gap, thin film side




740938 Quartz, 3.5 mm air gap, thin film bottom



740944 Quartz, thin film bottom




740945 Quartz, thin film side



Sample Tails


For Model 7404/7407

740932 Quartz to BN



shown with the 730937 Disposable BN cup


740935 Fiberglass



shown with the 730931 Kel-F bulk/powder upper and bottom cup


For Model 7410

740939 Quartz to Kel-F



shown with the 730931 Kel-F bulk/powder upper and bottom cup

740942 Quartz to BN




shown with the 730937 Disposable BN cup




730931 Kel-F bulk/powder upper and bottom cup




730933 Kel-F thin film side



730934 Kel-F thin film bottom




730935 Kel-F liquid upper and bottom cup



730937 Disposable BN cup



730938 BN thin film side



730939 BN thin film bottom



Recommended temperature usage

	Room temperature	Model 74034 oven	Model 74018 cryostat	Model 74035 variable temperature
Quartz	✓	✓	✓	✓
Fiberglass	✓			
Kel-F	✓		✓	✓ at RT and below
Boron nitride		✓		✓ above RT



Shipping Dimensions and Weight (w × d × h)

	Model 7404	Model 7407	Model 7410
Instrument console, electronics, head, and computer	122 cm × 84 cm × 165 cm (48 in × 33 in × 65 in) 392 kg (864 lb)	122 cm × 84 cm × 165 cm (48 in × 33 in × 65 in) 318 kg (700 lb)	122 cm × 84 cm × 165 cm (48 in × 33 in × 65 in) 318 kg (700 lb)
Electromagnet	109 cm × 94 cm × 135 cm (43 in × 37 in × 53 in) 471 kg (1038 lb)	122 cm × 97 cm × 128 cm (48 in × 38 in × 50 in) 860 kg (1896 lb)	107 cm × 107 cm × 114 cm (42 in × 42 in × 45 in) 1647 kg (3630 lb)
Electromagnet base	(magnet, base, and frame together)	(magnet, base, and frame together)	112 cm × 112 cm × 41 cm (44 in × 44 in × 16 in) 165 kg (363 lb)
Frame			122 cm × 109 cm × 147 cm (48 in × 43 in × 58 in) 341 kg (750 lb)
Power supply	(Included in instrument console)	109 cm × 79 cm × 117 cm (43 in × 31 in × 46 in) 331 kg (730 lb)	109 cm × 79 cm × 117 cm (43 in × 31 in × 46 in) 420 kg (926 lb)

Installation Dimensions and Weight (w × d × h)

	Model 7404	Model 7407	Model 7410
Instrument console, electronics, head, and computer	79 cm × 77 cm × 160 cm (31 in × 30 in × 63 in) 131 kg (289 lb)	79 cm × 77 cm × 160 cm (31 in × 30 in × 63 in) 57 kg (126 lb)	79 cm × 77 cm × 160 cm (31 in × 30 in × 63 in) 57 kg (126 lb)
Electromagnet, electromagnet base, and frame	84 cm × 82 cm × 140 cm (33 in × 32 in × 55 in) 307 kg (677 lb)	120 cm × 82 cm × 140 cm (47 in × 32 in × 55 in) 739 kg (1629 lb)	120 cm × 82 cm × 140 cm (47 in × 32 in × 55 in) 1392 kg (4259 lb)
Power supply	(Included in instrument console)	61 cm × 92 cm × 137 cm (24 in × 36 in × 54 in) 250 kg (551 lb)	61 cm × 92 cm × 137 cm (24 in × 36 in × 54 in) 273 kg (602 lb)

Application Notes

- Magnetic In-line Metrology for GMR Spin-Valve Sensors
- Finite Sample Size Effects on the Calibration of Vibrating Sample Magnetometers
- Low Moment Measurements with a Vibrating Sample Magnetometer
- Magnetic Anisotropy: Measurements with a Vector Vibrating Sample Magnetometer
- Measurements with a VSM—Permanent Magnet Materials
- The Performance of the Model 7400 VSM: Sensitivity
- Magnetic Media Measurements with a VSM

Visit www.lakeshore.com for the most up-to-date information

Site Requirements

A system-specific site prep checklist will be provided

Power

Instrumentation, computer, and optional vacuum pump require two standard single-phase electrical outlets (20 A maximum). Magnet power supply and optional recirculation chiller require 3-phase electrical outlets (21 A maximum).

Water

Electromagnet requires one supply and one return line for cooling with up to 15 L/min and 30 to 50 psi. Magnet power supply requires a minimum of 7.6 L/min with a maximum pressure of 80 psi and +15 °C to +30 °C water temperature.

Floor

The floor must support the weight of the magnet, supply, and the equipment used to move them into place. The weight of the console is negligible in comparison. Heavy concrete ground floors usually prove best, not only because they have the required strength, but such a floor also transmits minimal building vibration to the magnetometer.

The system also requires minimum spacing between each of the above three pieces and 0.75 m for access to the rear of the equipment. (See Installation Dimensions and Weight table).

Environment

The VSM requires a temperature-controlled environment that is relatively free of airborne dust and debris. There should be no equipment placed next to the VSM system that would emit or be susceptible to high levels of magnetic interference (distribution boxes, vibration equipment, x-ray machines, etc.)



Ordering Information

7400 Series Systems

7404	High sensitivity VSM with 4 in electromagnet, 643 magnet power supply
7407	High sensitivity VSM with 7 in electromagnet, 648 magnet power supply
7410	High sensitivity VSM with 10 in electromagnet, 648 magnet power supply

7400 Series Accessories

7404 and 7407 only	
740927	Sample tail kit, fiberglass to KelF [®] , RT, includes 1 each of 740933/4/5 and 730935, and 3 each of 730931/3/4
740928	Sample tail kit, quartz to BN, RT and oven, includes 1 each of 740930/1/2 and 730937
740929	Sample tail kit, fiberglass to KelF [®] and quartz, RT and cryogenic, includes 1 each of 740930/1
740930	One piece quartz sample tail/holder, RT and oven, thin film bottom
740931	One piece quartz sample tail/holder, RT and oven, thin film side
740932	Sample tail only, quartz to BN, oven, used with 730937/8/9 sample holder
740933	3.5 mm air gap, 1-piece fiberglass sample tail/holder, thin film side
740934	3.5 mm air gap, 1-piece fiberglass sample tail/holder, thin film bottom
740935	Sample tail only, fiberglass, used with 730931/3/4/5 sample holder
MMT-6J04-VG-06	4 in Hall probe for EM4, EM7 with 1 and 2 in coils (for 736 controller); replaces 735952 and 735954

7410 only

740939	Sample tail only, quartz to Kel-F [®] , RT, used with 730931/3/4 sample holder
740940	Sample tail kit, quartz to Kel-F [®] , RT, includes 1 each of 740937/8/9 and 730935, and 3 each of 730931/3/4
740941	Sample tail kit, quartz tail to BN cup, oven, includes 1 each of 740942 and 730937
740942	Sample tail only, quartz to BN, oven, used with 730937 sample holder
740943	Sample tail kit, fiberglass to Kel-F [®] and quartz, RT and cryogenic, includes 1 each of 740944/5
740944	1-piece quartz sample tail/holder, RT and oven, thin film bottom
740945	1-piece quartz sample tail/holder, RT and oven, thin film side
740937	3.5 mm air gap, 1-piece quartz sample tail/holder, thin film side
740938	3.5 mm air gap, 1-piece quartz sample tail/holder, thin film bottom
MMT-6J08-VG-06	8 in Hall probe for EM10 (for 736 controller); replaces 735958

7404, 7407, and 7410

730931	Sample holder cup, upper and bottom portion, Kel-F [®]
730933	Sample holder, thin film side, Kel-F [®]
730934	Sample holder, thin film bottom, Kel-F [®]
730935	Sample holder, liquid, upper and bottom portion, Kel-F [®]
730937	Sample holder, disposable, oven, BN
730938	Sample holder, thin film side, oven, BN
730939	Sample holder, thin film bottom, oven, BN
730904	Ceramic putty for oven sample mount
730907	Test sample sphere, NIST-traceable
730908	Test sample, 99% pure nickel sphere
730909	Test sample, 99% pure nickel 1 mm sphere

7400 Series Upgrades

Consult Lake Shore for information on upgrading your VSM system to the latest 7400 series

7400 Series Options

74018	Cryostat, variable temperature LN ₂ and LHe
74032	Vector option, Y coil for 2 in gap (requires large diameter pole caps with 740ESC X coils ordered separately)
74032XY	Vector option, X and Y coil set for 2 in gap (required large diameter pole caps ordered separately)
74033	Autorotation option, 736 controller (not field upgradable — VSM head must be returned to the factory)
74034	Oven, high temperature
74035	Single stage variable temperature
74046	Magnetoresistance probe — can be used with high and low temperature options (741-VTA not included)
741-VTA	Variable temperature option kit
730ESC	2 in pick-up coils
740EMSC	1 in standard pick-up coils

Magnetic Accessories

TPS-FRG-100/120V	Compact turbo pumping system; includes V-81 turbo pump (NW 40) with oil free dry scroll backing pump, FRG-700 full range gauge, controller, and interface cable to USB port; includes Agilent 24 month warranty NOTE: requires SYS-TP-KIT
TPS-FRG-220/240V-CE	Compact turbo pumping system; includes V-81 turbo pump (NW 40) with oil free dry scroll backing pump, FRG-700 full range gauge, controller, and interface cable to USB port; includes Agilent 24 month warranty NOTE: requires SYS-TP-KIT
SYS-TP-KIT	Includes all components necessary to connect NW 40 turbo pumping system to the vacuum port of any Lake Shore system (except probe stations)
E2M-110/120V	Two-stage rotary vacuum pump with mist filter; 110 to 120 VAC NOTE: requires SYS-RP-KIT
E2M-220/240V	Two-stage rotary vacuum pump with mist filter; 220 to 240 VAC NOTE: requires SYS-RP-KIT
SYS-RP-KIT	Includes all components necessary to connect E2M rotary pump to the vacuum port of any Lake Shore system (except probe stations)
1220-50	50 L LN ₂ Dewar with ½ in top withdraw port and 10 psi pressure relief valve
1230-60	60 L LHe Dewar

Got o www.lakeshore.com for the current list of available recirculating chillers

Lake Shore 7400 Series Vibrating Sample Magnetometers and Accessories

7404 VSM

7407 VSM

7410 VSM

Sample Rods and Tails

Options and Accessories

09/23/2015

Lake Shore Cryotronics, Inc.

575 McCorkle Blvd.

Westerville, OH 43082

Tel 614-891-2244

Fax 614-818-1600

info@lakeshore.com

www.lakeshore.com

